

PE Fast Chest CTA

Description:

Intended for PE studies and general Chest CTA protocols. Primary imaging factors include Temporal and S

Scan Settings

Head

Body



Primary

Large Body, Large, 80mm

0.5

1.531



SFOV, Focal Spot, Collimation

Rot Time

Pitch

Large Body, Large, 80mm

0.8

1.531



SFOV, Focal Spot, Collimation

Rot Time

Pitch

Large Body, Large, 80mm

0.8

0.992



SFOV, Focal Spot, Collimation

Rot Time

Pitch

Large Body, Large, 80mm

1.0

0.992



SFOV, Focal Spot, Collimation

Rot Time

Pitch

Large Body, XLarge, 80mm

1.0

0.508



SFOV, Focal Spot, Collimation

Rot Time

Pitch